

Archit PPH Cup: user guidelines

1. Kit includes plastic box, Archit PPH cup with two stoppers attached at two ports, Suction tube with connector, brush for cleaning purpose, small plastic container to keep formalin tablets and to be covered with gauze piece.
2. For insertion, **before/ after** C section on operation table or after normal delivery (After ruling out Vaginal wall trauma/ Cervical tear). Remove one stopper from one port. Fold the cup- "C" fold" or "7" fold. Apply lubricant jelly or lignocaine jelly. Insert the cup into lower/ mid vaginal canal (No need to push the high up). Release the fold. With one finger feel for the rim of the cup is opened. Apply the removed stopper.
3. If inserted **before** C section then keep the suction tube inserted from one port, attached to another suction machine. At another port stopper shall be in place. When needed intermittently suction is done as and when required. This helps in calculating the intraoperative blood loss accurately, keeps the operating field clean without much need to suction from uterine cavity. Prevents blood clots accumulating in the uterine cavity. It prevents soiling of OT table.
4. After shifting on bed. Remove one stopper and connect urine bag. Use another port for intermittent suction. Every 15-30 min for first 1-2 hours and later less frequently. Always do suction before removing the cup. Proper suction helps in removing the blood clots from the uterine cavity.
5. As per hospital protocol the PPH cup can be in situ till patient is mobilized out of bed.
6. Before removing the cup, Remove the urine catheter first. Talk to the patient, use any lubricant/ lignocaine jelly. After suction. Disconnect the urine bag. Keep one port open. Keep one finger inserted into lower vaginal canal. Gently pulldown the cup by side to side, up/ down movement. Once you feel broad funnel or rim, pinch it. This pinching will break the watertight suction between the rim and vaginal wall. Keep on gentle pull down words as above. And remove the PPH cup.
7. After use clean properly with any cleaning solution and brush as done with any other surgical instrument. Make it dry. Clean with spirit. Keep in the formalin box and it can be used for next patient.
8. Alternately any other disinfectant / sterilization method like Cidex solution or can be autoclaved. (Note: Frequent autoclaving cycles will cause hardening of the silicone material)

Note: This is the final patent granted product by the patent office of Government of India. All intellectual property rights are protected. Any unauthorized copy/ use is punishable as per law.

You can visit <https://architpphcup.com/video-intro/> for more information.

Any feedback is welcomed.